

# Systemic Rheumatic Disease

Identification of the antinuclear antibody (ANA) pattern remains a crucial step in the process of diagnosing the systemic rheumatic diseases. ANA patterns often give the clinician insight into which autoantibodies are present and indications of disease likelihood.

It has now been over 55 years since the LE cell was first described by Hargraves (1). During this time, substantial improvements in physician awareness, diagnostics and treatments have increased the survival rate for lupus patients from a 5 year rate of only 50% to a 10 year rate of >90% (2).

There is growing evidence that the appearance of autoantibodies may precede the onset of the disease, often by many years. Early detection of these autoantibodies may offer the opportunity for earlier diagnosis and treatment, improving the length and quality of life for the patient.

A recent study in the New England Journal of Medicine reported that, on average, at least one antinuclear antibody was present 3 years prior to the clinical diagnosis of SLE being made and at the time of diagnosis there were often 3 autoantibodies present (3). The earliest antibodies detected were ANA's, anti-SSA/Ro and anti-SSB/La antibodies. These antibodies were readily detected and reported using the HEp-2000® slide based assay.


The HEp-2000® substrate is Immuno Concepts' patented ANA substrate that has been consistently proven in independent studies to be superior to standard HEp-2 for the detection and identification of ANA's (4-7). When the unique SSA/Ro pattern is present on the HEp-2000® substrate the laboratory can immediately report the presence of anti-SSA/Ro antibodies to the clinician, potentially accelerating the correct diagnosis for the patient.

The photographs and charts presented here are designed to give the laboratorian a comprehensive guide to pattern recognition, antigen specificity, and disease association. Some of the more recently described patterns do not have disease association but are included for educational value.


### References:

- Hargraves, M., Richmond H., et al., Presentation of two bone marrow components, the tart cell and the LE cell. *Mayo Clin Proc.* 1948;27:25-28.
- Alamanos, Y., Voulgari, P. V., et al., Survival and mortality rates of systemic lupus erythematosus patients in northwest Greece. Study of a 21-year incidence cohort. *Rheumatology.* 2003;42(9):1122-1123.
- Arbuckle, M. R., McClain, M. T. et al. Development of autoantibodies before the clinical onset of systemic lupus erythematosus. *New England Journal of Medicine.* 2003;349(16):1526-1533.
- Pollock W, Toh BH. Routine immunofluorescence detection of Ro/SS-A autoantibody using HEp-2 cells transfected with human 60 kDa Ro/SS-A. *J.Clin.Pathol.* 1999;52:684-687.
- Bossuyt X, Meurs L., et al., Screening for autoantibodies to SS-A/Ro by indirect immunofluorescence using HEp-2000® cells. *Ann Clin Biochem.* 2000;37:216-219.
- Fritzler MJ, Hanson C., et al., Specificity of autoantibodies to SS-A/Ro on a transfected and overexpressed human 60 kDa Ro autoantigen substrate. *J.Clin.Lab.Anal.* 2002;16:103-108.
- Bossuyt X, Frans J., et al. Detection of Anti-SSA Antibodies by Indirect Immunofluorescence. *Clin Chem.* 2004;50(12):2361-2369.

# Flow Chart ANA Pattern and Confirmatory Testing



## ANA Negative



Immuno Concepts has taken care that the information and recommendations contained herein are accurate and compatible with current standards. Nevertheless, it is difficult to ensure that all of the information given is entirely accurate in all circumstances. Immuno Concepts disclaims any liability, loss or damage incurred as a consequence, directly or indirectly, of the use and application of any of the contents of this chart.

# Suggested Reading

### Testing Guidelines

- Kavanaugh A, Tomar R, Reveille J, Solomon DH, Homburger HA. Guidelines for clinical use of the antinuclear antibody test and tests for specific autoantibodies to nuclear antigens. *Arch.Pathol.Lab.Med.* 2000;124:71-81.
- Kavanaugh AF, Solomon DH, Amer Coll Rheumatology Ad Hoc C. Guidelines for immunologic laboratory testing in the rheumatic diseases: Anti-DNA antibody tests. *Arthritis & Rheumatism-Arthritis Care & Research.* Oct 15 2002;47(5):546-555.
- Solomon DH, Kavanaugh AJ, Schur PH, Amer Coll Rheumatology Ad Hoc C. Evidence-based guidelines for the use of immunologic tests: Antinuclear antibody testing. *Arthritis & Rheumatism-Arthritis Care & Research.* Aug 15 2002;47(4):434-444.
- Tozzoli R, Bizzaro N, Tonutti E, et al. Guidelines for the laboratory use of autoantibody tests in the diagnosis and monitoring of autoimmune rheumatic diseases. *Am.J.Clin.Pathol.* 2002;117:316-324.
- Reveille JD, Solomon DH, Amer Coll Rheumatology Ad Hoc C. Evidence-based guidelines for the use of immunologic tests: Anticentromere, Scl-70, and nucleolar antibodies. *Arthritis & Rheumatism-Arthritis Care & Research.* Jun 15 2003;49(3):399-412.


### HEp-2000®

- Keech CL, McCluskey J, Gordon TP. Transfection and overexpression of the human 60-kDa Ro/SS-A autoantigen in HEp-2 cells. *Clin.Immunol.Immunopathol.* 1994;73:146-151.
- Fritzler MJ, Miller BJ. Detection of autoantibodies to SS-A/Ro by indirect immunofluorescence using a transfected and overexpressed human 60 kD Ro autoantigen in HEp-2 cells. *J.Clin.Lab.Anal.* 1995;9:218-224.
- Pollock W, Toh BH. Routine immunofluorescence detection of Ro/SS-A autoantibody using HEp-2 cells transfected with human 60 kDa Ro/SS-A. *J.Clin.Pathol.* 1999;52:684-687.
- Bossuyt X, Meurs L, Mewis A, Mariën G, Blanckaert N. Screening for autoantibodies to SS-A/Ro by indirect immunofluorescence using HEp-2000® cells. *Ann Clin Biochem.* 2000;37:216-219.
- Fritzler MJ, Hanson C, Miller J, Eystathioy T. Specificity of autoantibodies to SS-A/Ro on a transfected and overexpressed human 60 kDa Ro autoantigen substrate. *J.Clin.Lab.Anal.* 2002;16:103-108.
- Reiff A, Haubruck H, Amos MD. Evaluation of a recombinant antigen enzyme-linked immunosorbent assay (ELISA) in the diagnostics of antinuclear antibodies (ANA) in children with rheumatic disorders. *Clinical Rheumatology.* May 2002;21(2):103-107.
- Bossuyt X, Frans J, Hendrickx A, et al. Detection of Anti-SSA Antibodies by Indirect Immunofluorescence. *Clin Chem.* 10 7 2004;50(12):2361-2369.



# Fluorescent Patterns

Associated with Antinuclear Antibody Detection



Speckled and SSA/Ro

Homogeneous and SSA/Ro


Centromere and SSA/Ro


SSA/Ro





# Antigen Chart Associated with ANA Patterns


Pattern observed by indirect immuno fluorescence	Type of antibody	Disease in which antibodies seen	Characteristics of antigenic determinants	Tests used to confirm specific antibody
<b>Homogeneous or peripheral</b>	Double-stranded DNA (dsDNA)	Characteristic of SLE, lower levels in other rheumatic diseases	Double-stranded DNA	IF using <i>Critidia luciliae</i> , RIA, ELISA, HA, CF
	Deoxyribonucleoprotein, soluble form	SLE, Drug-induced LE	DNA-histone complex	LE cell prep, RIA, ID, HA
	Histone	Drug-induced LE, SLE, others	Different classes of histone	ELISA, RIA
<b>Unusual Homogeneous</b>	Nuclear Envelope or Nuclear Membrane	Lupoid hepatitis	Nuclear Lamins A, B, C	none
<b>Speckled</b>	Sm	Marker antibody for SLE	Core proteins of the U1, U2, U4, U5, U6, U7, U11, and U12 snRNPs	ID, ELISA, HA, CF
	U1-RNP	High levels in MCTD and SLE, low levels in other diseases	A, C, and 70kDa proteins complexed with U1 RNA	ID, ELISA, HA, CF
	SSA/Ro	High prevalence in Sjogren's syndrome sicca complex, lower prevalence in other rheumatic diseases	52 and 60 kDa proteins complexed with Y1-Y5 RNA's	Confirmatory with HEp-2000® staining pattern; rule out other ENAs with ID, ELISA, CIEP, IP, western blot
	SSB/La	High prevalence in Sjogren's syndrome sicca complex, lower prevalence in other rheumatic diseases	48 kDa protein complexed with RNA polymerase III	ID ELISA, CIEP, IP, western blot
	Scl-70	Marker antibody for Scleroderma	70 kDa protein of topoisomerase	ID, ELISA
	PCNA	Marker antibody for SLE	Auxiliary protein to DNA polymerase δ	ID, CIEP
	Matrix	Seen in some patients with evolving rheumatic disease	Heterogeneous nuclear RNP (hnRNP), others	Confirmed by staining pattern
<b>Unusual Speckled</b>	NSp I, sp-100, MND, PBC 95	Some association with Primary Biliary Cirrhosis	95-100 kDa protein	none
	NSp II, CENP F	unknown	367 kDa protein	none
	Midbody	Low percentage of patients with scleroderma	Unidentified proteins of midbody region	none
	p 80 (coilin)	unknown	80 kDa protein in the coiled body	none
<b>Centromere</b>	Centromere	Seen in 57 - 82% of patients with CREST variant of scleroderma	CENP A, CENP B, CENP C,	Confirmed by staining pattern
<b>Nucleolar</b>	Clumpy Nucleolar	Scleroderma	Fibrillar, others?	none
	Speckled Nucleolar	Scleroderma	RNA polymerase I, others?	none
	Smooth Nucleolar	Polymyositis/Scleroderma	PM-1 (PM/Scl), others?	none


**Negative 1**  
  
 Report as: ANA Negative  
 Suspected antigen specificity: None  
 Clinical significance: None  
 Follow-up testing: None


**dsDNA positive 6**  
  
 Pattern: Positive staining of the kinetoplast in the *Critidia luciliae*  
 Report as: nDNA positive  
 Suspected antigen specificity: nDNA or dsDNA  
 Clinical significance: Marker antibody seen in 30 - 70% of SLE patients


**Homogeneous 2**  
  
 Report as: Homogeneous  
 Suspected antigen specificity: nDNA, DNP, Histone, DNA binding proteins, others?  
 Clinical significance: nDNA positive: marker antibody seen in 60% of SLE patients  
 Follow-up testing: Confirm nDNA antibodies


**dsDNA negative 7**  
  
 Pattern: No apparent staining of the kinetoplast in the *Critidia luciliae*  
 Report as: nDNA negative  
 Suspected antigen specificity: None


**Nuclear membrane 3**  
  
 Report as: Nuclear Membrane  
 Suspected antigen specificity: Nuclear lamins, others?  
 Clinical significance: May be seen in patients with SLE, RA, autoimmune hepatitis  
 Follow-up testing: None required

**Chromosomal Coating 8**  
  
 Report as: Suspect Chromosomal Coating Antibody  
 Pattern: Stains only the chromosomal area of metaphase mitotic cells  
 Report as: ANA negative  
 Suspected antigen specificity: Unknown  
 Clinical significance: Unknown  
 Follow-up testing: None required


**Nucleolar (clumpy) 4**  
  
 Report as: Nucleolar  
 Suspected antigen specificity: Fibrillar, others?  
 Clinical significance: May be seen in patients with SLE, RA, autoimmune hepatitis  
 Follow-up testing: None required


**Nucleolar (smooth) 9**  
  
 Report as: Nucleolar  
 Suspected antigen specificity: Pm/Scl  
 Clinical significance: Seen in patients with SLE, RA, polymyositis/scleroderma overlap  
 Follow-up testing: None required


**Nucleolar (speckled) 5**  
  
 Report as: Nucleolar  
 Suspected antigen specificity: RNA Polymerase I, NOR 90, others?  
 Clinical significance: Diffuse scleroderma  
 Follow-up testing: None required

**Sci-70 10**  
  
 Report as: Homogeneous, speckled, and nucleolar  
 Suspected antigen specificity: Topoisomerase I  
 Clinical significance: Topoisomerase I positive: marker antibody seen in 15 - 20% of patients with scleroderma  
 Follow-up testing: Confirm by ENA testing

**Speckled 11**  
  
 Report as: Speckled  
 Suspected antigen specificity: Sm, U1-RNP, others?  
 Clinical significance: Sm positive: marker antibody seen in 4 - 40% of SLE patients. RNP positive: seen in high titers in patients with MCTD and SLE, low titers in other diseases  
 Follow-up testing: Confirm by ENA testing

**Coarse Speckled 12**  
  
 Report as: Speckled  
 Suspected antigen specificity: U1-RNP, others?  
 Clinical significance: Sm positive: marker antibody seen in 4 - 40% of SLE patients. RNP positive: seen in high titers in patients with MCTD and SLE, low titers in other diseases  
 Follow-up testing: Confirm by ENA testing

**Speckled (flat or fine) 13**  
  
 Report as: Speckled  
 Suspected antigen specificity: SSA/Ro, SSB/La, others?  
 Clinical significance: SSA/Ro and SSB/La seen in high percentage of patients with primary Sjögren's syndrome, 30 - 40% of patients with SLE  
 Follow-up testing: Confirm by ENA testing


**HEp-2000® 200x 14**  
  
 A distinct speckled and nucleolar pattern seen in 10 - 20% of the interphase nuclei. These are the hyperexpressing cells. The remaining 80-90% of the interphase cells may or may not demonstrate staining. The chromosome region of the metaphase mitotic cells is negative. Original total magnification was 200x.

**HEp-2000® 400x 15**  
  
 Report as: SSA/Ro pattern  
 Suspected antigen specificity: SSA/Ro  
 Clinical significance: SSA/Rp positive: seen in 60-70% of patients with primary Sjögren's syndrome, 30 - 40% of patients with SLE  
 Follow-up testing: This pattern is confirmatory for SSA/Ro antibodies. Suggest ENA testing to rule out the presence of antibodies to other ENAs

# Fluorescent Patterns

## Associated with Antinuclear Antibody Detection

### Mitosis Stages of Cell Division



**Interphase**  
 DNA synthesis and chromosome replication occurs. Centrioles replicate.

**Prophase**  
 Chromosomes condense. Centrioles migrate to opposite ends of the cell. Nuclear membrane fragments.

**Metaphase**  
 The double-stranded chromosomes align at the center of the cell and the spindle fibers (generated from the centrioles) attach to the centromeres. Because the nuclear membrane has dissolved, nonchromosome specific antigens redistribute to outer region of the cell.


**Anaphase**  
 Chromosomes move to opposite poles.


**Telophase**  
 Spindle dissolves and nonchromosome specific antigens become contained within the chromosomal matrix as the nuclear membrane reforms.


**Cytokinesis**  
 Cell division is complete.


Note: chromosome staining demonstrated utilizing DNA specific antibody.


Tech Support 800-251-5115 [www.immunoconcepts.com](http://www.immunoconcepts.com)


**Nuclear Matrix 18**  
  
 Report as: Speckled  
 Suspected antigen specificity: hnRNP, others?  
 Clinical significance: Seen in patients with evolving rheumatic diseases, also Chronic Fatigue Syndrome  
 Follow-up testing: None required


**Midbody 23**  
  
 Report as: Atypical speckled, midbody  
 Suspected antigen specificity: Midbody  
 Clinical significance: May be seen in a low percentage of patients with scleroderma  
 Follow-up testing: None required, confirmed by staining pattern


**Cell Cycle Speckled 19**  
  
 Unusual Cell Cycle Dependent Speckled  
 Report as: Speckled or Unusual Speckled or Cell Cycle Dependent Speckled  
 Suspected antigen specificity: unknown  
 Clinical significance: Unknown  
 Follow-up testing: None required


**NspII 24**  
  
 NSp II (CENP F)  
 Report as: Atypical Speckled  
 Suspected antigen specificity: Centromere protein "F", others?  
 Clinical significance: Unknown, some association with malignancies  
 Follow-up testing: None required


**Centromere 20**  
  
 Report as: Centromere  
 Suspected Antigen Specificity: Centromere protein A, B, or C  
 Clinical significance: Seen in 57 - 96% of patients with the limited form of scleroderma (CREST)  
 Follow-up testing: None required, pattern confirmed by staining pattern


**Mixed ANA pattern 25**  
  
 Report as: Homogeneous and Speckled  
 Mixed ANA patterns can be caused by autoantibodies to several different antigens or by antibodies to an antigen located in several different areas of the cell nucleus. Titering and appropriate confirmatory testing for each pattern is recommended.


**Multiple Nuclear Dots 21**  
  
 Multiple Nuclear Dots (INSp I, PBC 95)  
 Report as: Atypical speckled  
 Suspected antigen specificity: 95 - 100 kD protein  
 Clinical significance: seen 27-44% of patients with PBC  
 Follow-up testing: None required


**Mixed ANA pattern 26**  
  
 Mixed ANA pattern on HEp-2000®  
 Report as: Homogeneous and SSA/Ro  
 Suspected antigen specificity: SSA/Ro—confirm by pattern; Homogeneous—dsDNA, histone or others?


**Multiple Nuclear Dots 21**  
  
 Multiple Nuclear Dots (INSp I, PBC 95)  
 Report as: Atypical speckled  
 Suspected antigen specificity: 95 - 100 kD protein  
 Clinical significance: seen 27-44% of patients with PBC  
 Follow-up testing: None required


**p80 coilin 22**  
  
 Report as: Atypical Speckled  
 Suspected antigen specificity: Coilin  
 Clinical significance: Unknown  
 Follow-up testing: None required


**Mixed ANA pattern 28**  
  
 Mixed ANA pattern on HEp-2000®  
 Report as: Centromere and SSA/Ro  
 Suspected antigen specificity: SSA/Ro—confirm by pattern; Centromere—confirm by pattern  
 Clinical significance: See significance for each respective pattern above


**suspect Mitochondrial 29**  
  
 Report as: ANA Negative, suspect Mitochondrial  
 Suspected antigen specificity: Mitochondria  
 Clinical significance: Seen in up to 90% of patients with PBC  
 Follow-up testing: Confirm on AMA specific substrate such as rodent tissue


**suspect Mitochondrial 29**  
  
 Report as: ANA Negative, suspect Mitochondrial  
 Suspected antigen specificity: Mitochondria  
 Clinical significance: Seen in up to 90% of patients with PBC  
 Follow-up testing: Confirm on AMA specific substrate such as rodent tissue


**Mitotic Spindle 34**  
  
 Report as: ANA Negative, Mitotic Spindle  
 Suspected antigen specificity: Mitotic spindles  
 Clinical significance: Seen in evolving rheumatic diseases  
 Follow-up testing: None required


**suspect Ribosomal P 30**  
  
 Report as: suspect Ribosomal P  
 Suspected antigen specificity: Ribosomal P proteins P0, P1, and P2  
 Clinical significance: Seen in up to 20% of patients with SLE  
 Follow-up testing: Confirm with anti-ribosomal specific assay


**NuMA 35**  
  
 Report as: ANA positive, Speckled, Mitotic spindle also present  
 Suspected antigen specificity: Nuclear Mitotic Apparatus  
 Clinical significance: Sjögren's Syndrome, PBC, others?  
 Follow-up testing: None required


**suspect Jo-1 31**  
  
 Report as: ANA Negative, cytoplasmic speckling  
 Suspected antigen specificity: Jo-1 (Histidyl tRNA synthetase)  
 Clinical significance: Seen in up to 35% of patients with myositis  
 Follow-up testing: Confirm with Jo-1 specific testing


**suspect Centriole 36**  
  
 Report as: ANA Negative, suspect Centriole  
 Suspected antigen specificity: Centriole proteins  
 Clinical significance: Unknown  
 Follow-up testing: None required

**suspect Endosome-GWB 32**  
  
 Report as: ANA Negative, cytoplasmic speckling  
 Suspected antigen specificity: Unknown  
 Clinical significance: Unknown  
 Follow-up testing: None required


**suspect Cytoskeletal 37**  
  
 Report as: ANA Negative, suspect Cytoskeletal  
 Suspected Antigen Specificity: Actin  
 Clinical significance: Seen in patients with autoimmune liver disease  
 Follow-up testing: Confirm on antigen specific substrate such as rodent tissue

**suspect Cytoskeletal 38**  
  
 Report as: ANA Negative, suspect Cytoskeletal  
 Suspected Antigen Specificity: tubulin, vimentin, others?  
 Clinical significance: Seen in patients with autoimmune liver disease  
 Follow-up testing: Confirm on antigen specific substrate such as rodent tissue

**Mechanical Damage 39**  
  
 Damage to the substrate caused by moving the coverslip or by touching the substrate with a pipet tip.

**Non-specific Film 40**  
  
 Non-specific film caused by adsorption of serum proteins to the slide surface.

**Heterophile Staining 41**  
  
 Report as: ANA Negative, cytoplasmic speckling  
 A result of non-specific antibody binding to substrate

**Evan's Blue counterstain 42**  
  
 Image demonstrates the effect of Evan's Blue counterstain