



## References

- Hargraves, M., Richmond H., et al., Presentation of two bone marrow components, the tart cell and the LE cell. *Mayo Clin Proc.* 1948;27:25-28.
- Alamanos, Y., Voulgari, P. V., et al., Survival and mortality rates of systemic lupus erythematosus patients in northwest Greece. Study of a 21-year incidence cohort. *Rheumatology.* 2003;42(9):1122-1123.
- Arbuckle, M. R., McClain, M. T. et al. Development of autoantibodies before the clinical onset of systemic lupus erythematosus. *New England Journal of Medicine.* 2003;349(16):1526-1533.
- Pollock W, Toh BH. Routine immunofluorescence detection of Ro/SS-A autoantibody using HEp-2 cells transfected with human 60 kDa Ro/SS-A. *J.Clin.Pathol.* 1999;52:684-687.
- Bossuyt X, Meurs L., et al., Screening for autoantibodies to SS-A/Ro by indirect immunofluorescence using HEp-2000® cells. *Ann Clin Biochem.* 2000;37:216-219.
- Fritzler MJ, Hanson C., et al., Specificity of autoantibodies to SS-A/Ro on a transfected and overexpressed human 60 kDa Ro autoantigen substrate. *J.Clin.Lab.Anal.* 2002;16:103-108.
- Bossuyt X, Frans J., et al. Detection of Anti-SSA Antibodies by Indirect Immunofluorescence. *Clin Chem.* 2004;50(12):2361-2369.



## Fluorescent Patterns

### Associated with Antinuclear Antibody Detection




## Suggested Reading

### Testing Guidelines

- Kavanaugh A, Tomar R, Reveille J, Solomon DH, Homburger HA. Guidelines for clinical use of the antinuclear antibody test and tests for specific autoantibodies to nuclear antigens. *Arch.Pathol.Lab.Med.* 2000;124:71-81.
- Kavanaugh AF, Solomon DH, Amer Coll Rheumatology Ad Hoc C. Guidelines for immunologic laboratory testing in the rheumatic diseases: Anti-DNA antibody tests. *Arthritis & Rheumatism-Arthritis Care & Research.* Oct 15 2002;47(5):546-555.
- Solomon DH, Kavanaugh AJ, Schur PH, Amer Coll Rheumatology Ad Hoc C. Evidence-based guidelines for the use of immunologic tests: Antinuclear antibody testing. *Arthritis & Rheumatism-Arthritis Care & Research.* Aug 15 2002;47(4):434-444.
- Tozzoli R, Bizzaro N, Tonutti E, et al. Guidelines for the laboratory use of autoantibody tests in the diagnosis and monitoring of autoimmune rheumatic diseases. *Am.J.Clin.Pathol.* 2002;117:316-324.
- Reveille JD, Solomon DH, Amer Coll Rheumatology Ad Hoc C. Evidence-based guidelines for the use of immunologic tests: Anticentromere, Scl-70, and nucleolar antibodies. *Arthritis & Rheumatism-Arthritis Care & Research.* Jun 15 2003;49(3):399-412.

### HEp-2000®


- Keech CL, McCluskey J, Gordon TP. Transfection and overexpression of the human 60-kDa Ro/SS-A autoantigen in HEp-2 cells. *Clin.Immunol.Immunopathol.* 1994;73:146-151.
- Fritzler MJ, Miller BJ. Detection of autoantibodies to SS-A/Ro by indirect immunofluorescence using a transfected and overexpressed human 60 kD Ro autoantigen in HEp-2 cells. *J.Clin.Lab.Anal.* 1995;9:218-224.
- Pollock W, Toh BH. Routine immunofluorescence detection of Ro/SS-A autoantibody using HEp-2 cells transfected with human 60 kDa Ro/SS-A. *J.Clin.Pathol.* 1999;52:684-687.
- Bossuyt X, Meurs L, Mewis A, Mariën G, Blanckaert N. Screening for autoantibodies to SS-A/Ro by indirect immunofluorescence using HEp-2000® cells. *Ann Clin Biochem.* 2000;37:216-219.
- Fritzler MJ, Hanson C, Miller J, Eystathioy T. Specificity of autoantibodies to SS-A/Ro on a transfected and overexpressed human 60 kDa Ro autoantigen substrate. *J.Clin.Lab.Anal.* 2002;16:103-108.
- Reiff A, Haubruck H, Amos MD. Evaluation of a recombinant antigen enzyme-linked immunosorbent assay (ELISA) in the diagnostics of antinuclear antibodies (ANA) in children with rheumatic disorders. *Clinical Rheumatology.* May 2002;21(2):103-107.
- Bossuyt X, Frans J, Hendrickx A, et al. Detection of Anti-SSA Antibodies by Indirect Immunofluorescence. *Clin Chem.* 2004;50(12):2361-2369.



Speckled and SSA/Ro



Homogeneous and SSA/Ro



Centromere and SSA/Ro



SSA/Ro

## Systemic Rheumatic Disease

Identification of the antinuclear antibody (ANA) pattern remains a crucial step in the process of diagnosing the systemic rheumatic diseases. ANA patterns often give the clinician insight into which autoantibodies are present and indications of disease likelihood.


It has now been over 55 years since the LE cell was first described by Hargraves (1). During this time, substantial improvements in physician awareness, diagnostics and treatments have increased the survival rate for lupus patients from a 5 year rate of only 50% to a 10 year rate of >90% (2).

There is growing evidence that the appearance of autoantibodies may precede the onset of the disease, often by many years. Early detection of these autoantibodies may offer the opportunity for earlier diagnosis and treatment, improving the length and quality of life for the patient.

A recent study in the New England Journal of Medicine reported that, on average, at least one antinuclear antibody was present 3 years prior to the clinical diagnosis of SLE being made and at the time of diagnosis there were often 3 autoantibodies present (3). The earliest antibodies detected were ANA's, anti-SSA/Ro and anti-SSB/La antibodies. These antibodies are readily detected and reported using the HEp-2000® slide based assay.

The HEp-2000® substrate is Immuno Concepts' patented ANA substrate that has been consistently proven in independent studies to be superior to standard HEp-2 for the detection and identification of ANA's (4-7). When the unique SSA/Ro pattern is present on the HEp-2000® substrate the laboratory can immediately report the presence of anti-SSA/Ro antibodies to the clinician, potentially accelerating the correct diagnosis for the patient.

### ANA Positive




### ANA Negative

If the ANA is negative, but disease is suspected, test for specific autoantibodies

## Flow Chart

### ANA Pattern and Confirmatory Testing



Immuno Concepts has taken care that the information and recommendations contained herein are accurate and compatible with current standards. Nevertheless, it is difficult to ensure that all of the information given is entirely accurate in all circumstances. Immuno Concepts disclaims any liability, loss or damage incurred as a consequence, directly or indirectly, of the use and application of any of the contents of this chart.

# Antigen Chart


## Associated with ANA Patterns

Pattern observed by indirect immuno fluorescence	Type of antibody	Disease in which antibodies seen	Characteristics of antigenic determinants	Tests used to confirm specific antibody
<b>Homogeneous or peripheral</b>	Double-stranded DNA (nDNA, dsDNA)	Characteristic of SLE, lower levels in other rheumatic diseases	Double-stranded DNA	IF using Crithidia luciliae, RIA, ELISA, HA, CF
	Deoxynucleoprotein, soluble form	SLE, Drug-induced LE	DNA-histone complex	LE cell prep, RIA, ID, HA
	Histone	Drug-induced LE, SLE, others	Different classes of histone	ELISA, RIA
<b>Unusual Homogeneous</b>	Nuclear Envelope or Nuclear Membrane	Lupoid hepatitis	Nuclear Lamins A, B, C	none
<b>Speckled</b>	Sm	Marker antibody for SLE	Core proteins of the U1, U2, U4, U5, U6, U7, U11, and U12 snRNPs	ID, ELISA, HA, CF
	U1-RNP	High levels in MCTD and SLE, low levels in other diseases	A, C, and 70kDa proteins complexed with U1 RNA	ID, ELISA, HA, CF
	SSA/Ro	High prevalence in Sjogren's syndrome sicca complex, lower prevalence in other rheumatic diseases	52 and 60 kDa proteins complexed with Y1-Y5 RNA's	Confirmatory with HEp-2000® staining pattern; rule out other ENAs with ID, ELISA, CIEP, IP, western blot
	SSB/La	High prevalence in Sjogren's syndrome sicca complex, lower prevalence in other rheumatic diseases	48 kDa protein complexed with RNA polymerase III	ID ELISA, CIEP, IP, western blot
	Scl-70	Marker antibody for Scleroderma	70 kDa protein of topoisomerase	ID, ELISA
	PCNA	Marker antibody for SLE	Auxiliary protein to DNA polymerase δ	ID, CIEP
	Matrix	Seen in some patients with evolving rheumatic disease	Heterogeneous nuclear RNP (hn RNP), others	Confirmed by staining pattern
<b>Unusual Speckled</b>	Nsp I, sp-100, MND, PBC 95	Some association with Primary Biliary Cirrhosis	95-100 kDa protein	none
	Nsp II, CENP F	unknown	367 kDa protein	none
	Midbody	Low percentage of patients with scleroderma	Unidentified proteins of midbody region	none
	p 80 (coilin)	unknown	80 kDa protein in the coiled body	none
<b>Centromere</b>	Centromere	Seen in 57 - 82% of patients with CREST variant of scleroderma	CENP A, CENP B, CENP C,	Confirmed by staining pattern
<b>Nucleolar</b>	Clumpy Nucleolar	Scleroderma	Fibrillarlin, others?	none
	Speckled Nucleolar	Scleroderma	RNA polymerase I, others?	none
	Smooth Nucleolar	Polymyositis/Scleroderma	PM-1 (PM/Scl), others?	none




Negative 1

Report as: ANA Negative  
 Suspected antigen specificity: None  
 Clinical significance: None  
 Follow-up testing: None



Speckled 2

Report as: Speckled  
 Suspected antigen specificity: Sm, U1-RNP, others?  
 Clinical significance: Sm positive: marker antibody seen in 4 - 40% of SLE patients. RNP positive: seen in high titers in patients with MCTD and SLE, low titers in other diseases.  
 Follow-up testing: Confirm by ENA testing



Homogeneous 3

Report as: Homogeneous  
 Suspected antigen specificity: nDNA, DNP, Histone, DNA binding proteins, others?  
 Clinical significance: nDNA positive: marker antibody seen in 60% of SLE patients  
 Follow-up testing: Confirm nDNA antibodies




Nucleolar (clumpy) 4

Report as: Nucleolar  
 Suspected antigen specificity: Fibrillarlin, others?  
 Clinical significance: May be seen in patients with systemic sclerosis.  
 Follow-up testing: None required




HEp-2000® 200x 5

A distinct speckled and nucleolar pattern seen in 10 - 20% of the interphase nuclei. These are the hyperexpressing nuclei. The remaining 80-90% of the interphase cells may or may not demonstrate staining. The chromosome region of the metaphase mitotic cells is negative. Original total magnification was 200x.



HEp-2000® 400x 6

Report as: SSA/Ro pattern  
 Suspected antigen specificity: SSA/Ro  
 Clinical significance: SSA/Ro positive: seen in 60-70% of patients with primary Sjögren's syndrome, 30 - 40% of patients with SLE.  
 Follow-up testing: This pattern is confirmatory for SSA/Ro antibodies. Suggest ENA testing to rule out the presence of antibodies to other ENAs



Centromere 7

Report as: Centromere  
 Suspected Antigen Specificity: Centromere protein A, B, or C  
 Clinical significance: Seen in 57 - 96% of patients with the limited form of scleroderma (CREST).  
 Follow-up testing: None required, pattern confirmed by staining pattern




Scl-70 8

Report as: Homogeneous, speckled, and nucleolar  
 Suspected antigen specificity: Topoisomerase I  
 Clinical significance: Topoisomerase I positive: marker antibody seen in 15 - 20% of patients with scleroderma.  
 Follow-up testing: Confirm by ENA testing




PCNA 9

Report as: Speckled, possible PCNA  
 Suspected antigen specificity: DNA polymerase δ (cyclin)  
 Clinical significance: PCNA positive: marker antibody seen in 2 - 10% of SLE patients  
 Follow-up testing: Confirm by ENA testing




Mixed ANA pattern 10

Report as: Homogeneous and Speckled  
 Mixed ANA patterns can be caused by autoantibodies to several different antigens or by antibodies to an antigen located in several different areas of the cell nucleus. Titration and appropriate confirmatory testing for each pattern is recommended.



Mixed ANA pattern 11

Mixed ANA pattern on HEp-2000®  
 Report as: Homogeneous and SSA/Ro  
 Suspected antigen specificity: SSA/Ro—confirm by pattern; Homogeneous—dsDNA, histone or others?



Mixed ANA pattern 12

Mixed ANA pattern on HEp-2000®  
 Report as: Speckled and SSA/Ro  
 Suspected antigen specificity: SSA/Ro—confirm by pattern; Speckled—Sm, U1-RNP, SSA/Ro, SSB/La, others