# IMMUNO CONCEPTS Optional Powder Components

**For In Vitro Diagnostic Use** IVD

**For Professional Use**

**INTENDED USE:** This component is for use in the Immuno Concepts Test Systems.

This component has been tested in the Immuno Concepts Test Systems and found to meet performance requirements; it has not been tested for performance in other manufacturer's test systems.

Use: PBS buffer and Colorzyme® color reagent must be dissolved in deionized or distilled water before use.

Storage: All components can be stored at 2-25 °C.

Stability: All components remain stable at least 12 months from date of manufacture. Do not use any component beyond its expiration date.

**PBS Buffer Powder: Kit Catalog No. 1012 (20)**

**Non-Reactive Component: PBS Buffer Powder, Catalog No. 1011.** Phosphate-buffered saline powder (0.01 M, pH 7.4 ± 0.2). Each pouch contains sufficient buffer powder to make 1 liter. (One pouch of buffer powder is supplied for each five slides in complete test kits.)

**Preparation:** Dissolve one pouch of buffer powder in 1 liter of deionized or distilled water, cover, and store between 2-25 °C for up to four weeks or until signs of contamination or other visible changes occur.

**Color Reagent Powder: Catalog No. 4066 (5)**

**Reactive Component: Color Reagent Powder, Catalog No. 4066.** HRP specific enzyme substrate powder, containing 4 chloro 1-naphthol. Each package contains powder to make 150 ml of self-activating Colorzyme® color reagent.

**Preparation:** Dissolve contents of one pouch in 150 ml of deionized or distilled water. Mix well until completely dissolved. Color reagent is stable for 30 days at room temperature in a closed container. Color reagent may be reused for up to 30 days or until any color change or precipitate is visible. Cloudiness or opalescence, with no visible precipitate upon reuse is normal. Depending on the rate of use, 150 ml of Colorzyme® color reagent can be used with up to twenty slides.

**FOR TECHNICAL ASSISTANCE: 1-800-251-5115 (Toll Free) or e-mail: [technicalsupport@immunoconcepts.com](mailto:technicalsupport@immunoconcepts.com)**

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**MDSS GmbH**