FITC QUALITY CONTROL SLIDE™

Optimal fluorescent microscopy requires consistent performance of the microscope optical systems. Results obtained with serum reference standards, such as those available from the Center for Disease Control, Atlanta, Georgia, may be useful indicators of changes in sensitivity, but they were designed more as guides to identify staining patterns. The performance characteristics of the fluorescent microscope itself are an important factor in the evaluation of fluorescent intensity and in the subsequent assessment of sensitivity. Comparison of fluorescent staining intensities on fluorescent microscopes requires the use of quantitative standards that have the same fluorescent spectra as the cells used in the FITC test procedures.

The FITC QC Slide™ permits daily monitoring of microscope performance and comparison of test results requiring fluorescent microscopy. In addition, these reference slides serve as a convenient and sensitive check on field illumination uniformity, as well as stability of the illuminator source. It is highly recommended that these slides be used in conjunction with specific antibody reference standards. The relative brightness of cells stained in FITC test procedures can be compared to these reference microbeads.

The Immuno Concepts FITC QC Slide™ contains five levels of fluorescent microbead standards. The immobilized microbeads contain uniform levels of fluorescein isothiocyanate (FITC) standardized to demonstrate the same excitation and emission spectra as fluorescence read from cells complexed with fluorescent antibody conjugates. The intensities stated were set in our laboratory in collaboration with Dr. Eng Tan at the Scripps Clinic and Research Foundation. The stated intensities may differ from yours. To be assured of a minimum threshold of sensitivity with your fluorescent microscope system, the endpoint of positive staining should be observed with the microbeads on the 1+ or +/- wells. If you do not obtain results in this range, see troubleshooting sensitivity differences below.

USE OF THE IMMUNO CONCEPTS FITC QUALITY CONTROL SLIDE™

1. Set up the microscope for viewing at a total magnification of 200X.
2. Begin viewing at the +/- well and proceed to the 1+, 2+, 3+ and 4+ well, noting at which well positive staining of the microbeads is observed. (Move quickly from field to field to reduce photo bleaching.) The endpoint is the first well in which positive staining in the majority of the microbeads is observed.
3. Store the slides at 2-10°C and protect from light when not in use.

NOTE: With repeated use of the Immuno Concepts FITC QC Slide™, areas of diminished fluorescence may be observed. Always view fields with consistent intensity between beads. Discard the slide when over 50% of the microbeads in each field demonstrate diminished fluorescence.

TROUBLESHOOTING SENSITIVITY DIFFERENCES

Low Sensitivity:

If your endpoint for positive staining of the fluorescent microbeads is 4+, 3+, or 2+, check the following items:

1. **Proper alignment of the microscope light path.** Refer to your fluorescent microscope manual for instructions.
2. **The type of light source used.** Mercury light sources will produce greater excitation energy at 495 nm than Quartz/Halogen. The 50-watt, 100-watt, and 200-watt Mercury light sources differ little in excitation energy at 495nm. The 100-watt Quartz/Halogen light sources will produce greater excitation energy at 495nm than 50-watt Quartz/Halogen.
3. **The condition and age of the light source.** This is particularly true for Mercury light sources which generally exhibit a gradual reduction in excitation energy at 495nm prior to burning out. This gradual reduction in excitation energy can result in a significant loss in sensitivity over several weeks. This problem can be avoided by keeping a time log. For best results replace 50-watt mercury bulbs at 100 hours, and 100 or 200-watt mercury bulbs at 200 hours. Quartz/Halogen light sources generally do not exhibit a gradual reduction in excitation energy prior to burning out.
4. **The type of exciter filter used.** Interference exciter filters provide greater sensitivity over a much narrower wavelength than absorption exciter filters. Refer to your fluorescent microscope manual or sales representative for more information.
5. **The numerical aperture of the objective.** With incident light fluorescence (Epi), fluorescence is increased exponentially as the numerical aperture (NA) of the objective is increased additively. This may cause a 40X objective with a NA of 0.65 to read one or more dilutions lower than a 40X objective with a NA of 0.85. The numerical aperture is printed on the side of the objective. Sensitivity of transmitted light fluorescent microscopy is not affected by NA.
6. **Reference sera and normal studies.** Sensitivity is always relative to performance characteristics of the test with positive reference sera and your specific patient population.

High Sensitivity:

If strong positive staining of the fluorescent microbeads is seen at the +/- well, check the following items:

1. **Reference sera and normal studies.** Sensitivity is always relative to performance characteristics of the test with positive reference sera and your specific patient population.
2. **Suppression filters to reduce sensitivity.** Suppression filters reduce specific wavelengths of excitation and may be useful in reducing sensitivity. Refer to your fluorescent microscope manual or sales representative for more information.

**For Technical Assistance: telephone 1-916-363-2649 or 1-800-251-5115 (USA Toll Free)
Facsimile: 1-916-363-2843 or email technicalsupport@immunoconcepts.com**

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