IMMUNO CONCEPTS ANA Optional Components

For in vitro diagnostic use **IVD**

Component, Composant, Componente, Komponente:

Cat#, REF

Control Sera / Sérums de Contrôle / Sieri di Controllo
Sueros de Control / Kontrollseren / Kontrollsera:

2021	2022	2023	2025	2026	2027
2029	2031	2035	2037	2041	

Enzyme antibody reagents / Réactif Immuno-Enzymatique
Reagente Anticorpo Enzimatico / Anticuerpo Conjugado con Enzima
Enzym-Antikörperreagenz / Enzymantintikroppreagens:

4075 4075-Ro

INTENDED USE: This component is for use in the Immuno Concepts ANA/ANA-Ro Test Systems.

This component has been tested in the Immuno Concepts ANA/ANA-Ro Test Systems and found to meet performance requirements; it has not been tested for performance in other manufacturer's test systems.

Use: All components come ready to use with no aliquoting or reconstitution required (except the Titratable Control Serum which is treated as an undiluted patient sample.).

Storage: All components can be stored under refrigeration at 2-10°C.

Stability: All components remain stable at least 12 months from date of manufacture. Do not use any component beyond its expiration date.

PRECAUTIONS

1. All human source materials used for this product have been tested and found to be negative (not repeatedly reactive) for antibodies to Human Immunodeficiency Virus-1 (HIV-1), Human Immunodeficiency Virus-2 (HIV-2), hepatitis C virus (HCV) and for hepatitis B surface antigen (HBsAg) by FDA approved methods. However, no test method can offer complete assurance that HIV-1, HIV-2, hepatitis C, hepatitis B, or other infectious agents are absent. Therefore, all kit materials should be handled in the same manner as potentially infectious material.

2. Sodium azide (0.09%) is used as a preservative. When disposing of reagents, flush with ample volumes of tap water to prevent potential residues in plumbing. Sodium azide is a poison and may be toxic if ingested.

FOR TECHNICAL ASSISTANCE: 1-800-251-5115 (Toll Free) or e-mail: technicalsupport@immunoconcepts.com

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MDSS GmbH