IMMUNO CONCEPTS AUTO I.D.® Optional Components

For in vitro diagnostic use **IVD**

INTENDED USE: This component is for use in the Immuno Concepts Auto ID® Test Systems.

This component has been tested in the Immuno Concepts Auto ID® Test Systems and found to meet performance requirements; it has not been tested for performance in other manufacturer's test systems.

USE: All control sera come ready to use with no diluting, aliquoting, or reconstitution required. The nuclear antigen comes lyophilized and must be reconstituted with distilled or deionized water prior to use.

STORAGE: All components can be stored in the refrigerator at 2-10°C. Reconstituted nuclear antigen should be used within 5 days or stored frozen in aliquots at -20°C or less. Do not freeze AUTO I.D.® plates

STABILITY: All control sera are stable for at least 12 months from date of manufacture. Agarose plates are stable for 24 months from date of manufacture. Lyophilized nuclear antigen is stable at least 12 months from date of manufacture. After reconstitution, nuclear antigen is stable for 5 days at 2-10°C, or 90 days frozen at -20°C or less. For best results store reconstituted antigen in 30 µl aliquots at -20°C or less. Avoid repeated freeze/thawing of reconstituted nuclear antigen.

Positive Control Serum (0.5 ml):

USE: All control sera come ready to use with no diluting, aliquoting, or reconstitution required.

Cat. No: Ready-to-use vial containing:

6001	Human antibodies reactive with U1-Ribonucleoprotein (RNP) nuclear antigen
6002	Human antibodies reactive with Smith (Sm) and U1-Ribonucleoprotein (RNP) nuclear antigens
6003	Human antibodies reactive with Smith (Sm) nuclear antigen
6004	Human antibodies reactive with Jo-1 antigen
6005	Human antibodies reactive with Scl-70 nuclear antigen
6006	Human antibodies reactive with proliferating cell nuclear antigen (PCNA)
7001	Human antibodies reactive with SSA/Ro nuclear antigen
7002	Human antibodies reactive with SSA/Ro and SSB/La nuclear antigens
7003	Human antibodies reactive with SSB/La nuclear antigen

AUTO I.D.® Nuclear Antigen Cat. No. 6050 (0.2 ml): Lyophilized mammalian extracted nuclear antigen containing Smith (Sm), U1-Ribonucleoprotein (RNP), SSA/Ro, SSB/La, and other antigens that commonly react with antibodies from patients with systemic rheumatic disease. Each vial must be reconstituted with 200 µl of distilled or deionized water prior to use. Always allow freshly reconstituted nuclear antigen to sit at least 30 minutes at room temperature prior to use to make sure antigen is completely dissolved.

Preparation: Remove the metal cap from the nuclear antigen vial. Carefully lift the rubber stopper to vent the vial. Remove the stopper and add 200 µl of distilled or deionized water to the vial. Replace the rubber stopper and gently swirl to dissolve contents. Allow reconstituted antigen to sit at least 30 minutes prior to use to insure antigen is completely dissolved. Reconstituted antigen may appear turbid or cloudy. Swirl immediately before use. For best results store reconstituted antigen in 30 µl aliquots at -20°C or less. Allow antigen to reach room temperature before use.

AUTO I.D.® Plates Cat. No. 7010: Seven-well agarose plates that include prenumbered wells to ensure easy identification of patients.

Preparation: Allow plate to reach room temperature (18-24°C) prior to opening foil pouch. Carefully remove plate from foil pouch. Condensation present on the inside lid of the plate can be removed with bibulous paper or lint-free paper toweling. Avoid touching the agarose.

PRECAUTIONS

1. All human source materials used for this product have been tested and found to be negative (not repeatedly reactive) for antibodies to Human Immunodeficiency Virus-1 (HIV-1), Human Immunodeficiency Virus-2 (HIV-2), hepatitis C virus (HCV) and for hepatitis B surface antigen (HBsAg) by FDA approved methods. However, no test method can offer complete assurance that HIV-1, HIV-2, hepatitis C, hepatitis B, or other infectious agents are absent. Therefore, all kit materials should be handled in the same manner as potentially infectious material.

2. Sodium azide (0.1%) is used as a preservative. When disposing of reagents, flush with ample volumes of tap water to prevent potential residues in plumbing. Sodium azide is a poison and may be toxic if ingested.

FOR TECHNICAL ASSISTANCE: 1-800-251-5115 (Toll Free) or e-mail: technicalsupport@immunoconcepts.com

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MDSS GmbH